

Step-by-step manual for APEX Pipeline in OT-2

Martyna Kasprzyk¹, Michael Herrera², and Giovanni Stracquadanio^{1,*}

¹School of Biological Sciences, The University of Edinburgh, Edinburgh EH9 3BF, United Kingdom

²School of Chemistry, The University of Edinburgh, Edinburgh EH9 3BF, United Kingdom

*Corresponding author. Email: giovanni.stracquadanio@ed.ac.uk.

Contents

Getting Started	2
Downloading and Running the Pipeline	4
Configuration and Setup	4
Protocol 1: Heat Shock Transformations	7
Protocol 2: Colony Selection	10
Protocol 3: Colony Sampling	12
Protocol 4: Protein Expression Induction	15
Protocol 5: Colony PCR Sample Preparation	18
Protocol 6: SDS PAGE Sample Preparation	21
Protocol 7: OD₆₀₀ Reading Sample Preparation	24

Getting Started

This manual provides all the necessary information to set up and operate the Automated Protein Expression (APEX) pipeline using the OT-2. In *Escherichia coli* (*E. coli*), a standard protein expression workflow typically consists of four sequential steps. We've designed APEX with a modular architecture that mirrors this process, implementing four self-contained modules, each describing a specific OT-2 protocol.

Protocol 1 Heat shock transformation to deliver plasmid DNA into competent cells.

Protocol 2 Antibiotic selection on agar plates to isolate colonies carrying the DNA payload.

Protocol 3 Colony sampling and culturing of successfully transformed cells.

Protocol 4 Chemical induction of protein expression.

Additionally, we have developed supplementary protocols for:

Protocol 5 Colony PCR sample preparation for rapid screening of successful transformants.

Protocol 6 SDS-PAGE sample preparation for protein expression analysis and verification.

Protocol 7 Optical Density at 600 nm (OD₆₀₀) measurement sample preparation for monitoring bacterial growth.

APEX operates through a structured configuration system requiring two key input files: a JavaScript Object Notation (JSON) file for OT-2 hardware settings and a Comma Separated Value (CSV) file containing experiment-specific liquid handling instructions. These files provide the essential parameters for laboratory protocol setup and execution.

The system integrates these inputs into a master Python script that can be directly imported into the robot via the Opentrons app. Simultaneously, APEX generates a PDF document featuring both textual instructions and visual representations of the labware layout, serving as a practical reference guide for users setting up experiments (Figure 1).

Opentrons Tough 96 Well Plate 200 µL PCR Full Skirt in slot thermocycle




Figure 1: Example of labware visualised by APEX.

For validation purposes, APEX includes simulates each protocol that produces a detailed execution log, similar to the run log in the Opentrons App. This simulation output allows to preview and verify all protocol steps before physical implementation.

The pipeline is managed using Nextflow, a workflow manager that streamlines the execution of complex computational processes. Furthermore, APEX is fully containerised using Docker, which encloses all software dependencies into a single package. This containerisation ensures consistent deployment across any computational environment without requiring additional software installations.

System Requirements

Installing Java

Ensure Java 8 or later is installed on your system. Visit <https://www.java.com/en/download/> for installation instructions if needed. Java is required to run Nextflow.

Installing Nextflow

Open a terminal and execute the following command to download and install Nextflow:

```
curl -s https://get.nextflow.io | bash
```

Move the Nextflow executable to a directory in your PATH for easier access:

```
sudo mv nextflow /usr/local/bin
```

Verify the installation:

```
nextflow -version
```

Downloading and Running the Pipeline

The APEX pipeline is available on GitHub. Use either of the following commands to download it:

```
nextflow pull stracquadaniolab/apex-nf
```

Or directly with the GitHub URL:

```
nextflow pull https://github.com/stracquadaniolab/apex-nf
```

The pipeline will be downloaded to:

```
$HOME/.nextflow/assets/stracquadaniolab/apex-nf
```

To run the complete pipeline:

```
nextflow run stracquadaniolab/apex-nf
```

To run a specific protocol (e.g., Protocol 1):

```
nextflow run stracquadaniolab/apex-nf -entry PROTOCOL_1
```

To view project details, installation location or check for updates:

```
nextflow info stracquadaniolab/apex-nf
```

Configuration and Setup

The pipeline is organised in the "assets" folder with the following structure:

- testdata** JSON and CSV example files for protocol configurations and experimental data.
- instructions** R markdown template files containing experiment setup instructions.
- labware** JSON files defining supported labware. Custom labware definitions can be added here (more information here).
- protocols** Python template files for each protocol.
- results** Output files including ready-to-use Python protocols and PDF instructions.

Input Files

After downloading the necessary software, setting up your experiment with APEX is straightforward. Create a JSON and a CSV file containing all the details for your experiment. These files are critical as they dictate the setup and execution of your protocol.

Example files can be found in the `testdata` folder of the APEX directory once it has been downloaded. These examples serve as templates for proper formatting and contain all necessary parameters. The input files follow these naming conventions:

- **JSON File:** Name your file as `protocol-*-config.json`, where `*` represents the protocol ID number. For example, a transformation experiment would use `protocol-1-config.json`.
- **CSV File:** Similarly, name your CSV file as `protocol-*-data.csv`. For transformation, this would be `protocol-1-data.csv`.

Understanding JSON Configuration Files

JavaScript Object Notation (JSON) configuration file is used in the pipeline for setting up the Open-
trons. This file includes the details about the hardware, such as types of pipettes available and layout of the Open-
trons deck, specifying where each piece of labware should be placed.

Here's an example JSON configuration file that outlines a simple setup for an OT-2 robot:

```
{
  "pipette_name": "p20_multi_gen2",
  "source_plate_name": "armadillo_96_wellplate_200ul_pcr_full_skirt",
  "source_plate_slot": 4
}
```

For each protocol within the pipeline, this manual provides a specific JSON template with detailed explanations of all parameters.

To create a JSON file, you can use for example a simple text editor like Notepad (Windows) or TextE-
dit (Mac, set to plain text mode). Type your JSON data and save the file with a `.json` extension by selecting Save As from the File menu. Move the saved file to the `testdata` folder within the APEX project directory.

Understanding CSV Data Files

Comma Separated Value (CSV) files provide detailed instructions for operations like liquid transfers within the protocol. They allow users to specify parameters such as transfer volumes, source and destination wells, and other experimental variables without modifying the underlying code.

Example of a CSV file:

```
source_well,destination_well,volume  
A1,A2,50  
B1,C2,25
```

For each protocol within the pipeline, this manual provides a specific CSV template with detailed explanations of all parameters.

To create a CSV file, use either a basic text editor or a spreadsheet program like Microsoft Excel. Enter your data, ensuring each field is separated by a comma without spaces. Save the file with a '.csv' extension, selecting 'CSV (Comma delimited) (*.csv)' if using a spreadsheet program. Ensure the file is saved in the `testdata` folder within the APEX project directory.

Output files

After running the pipeline, all outputs will be saved in the `results` folder. Output files follow these naming conventions:

- **Python protocol:** `protocol-*.py`, where `*` represents the ID number of the protocol. For example, the Python protocol for a transformation experiment would be named `protocol-1.py`.
- **Protocol simulation:** `protocol-*.txt`. For the transformation example, this would be `protocol-1.txt`.
- **PDF instructions:** `protocol-*.pdf`. For the transformation example, this would be `protocol-1.pdf`.

Protocol 1: Heat Shock Transformations

Design

Use the JSON (protocol-1-config.json) and CSV (protocol-1-data.csv) files to create a protocol tailored for heat shock transformations. Define parameters as described in Table 1 and Table 2 below.

Parameter	Description	Example
right_pipette_name	Model of the pipette on the right mount (specify "NA" if not used)	"p20_multi_gen2"
right_pipette_tiprack_name	Tip rack type compatible with right pipette	"opentrons_96_tiprack_20ul"
right_pipette_tiprack_slot	Deck slot(s) where right pipette tip racks are placed	[4, 5]
left_pipette_name	Model of the pipette on the left mount (specify "NA" if not used)	"p300_multi_gen2"
left_pipette_tiprack_name	Tip rack type compatible with left pipette	"opentrons_96_tiprack_300ul"
left_pipette_tiprack_slot	Deck slot(s) where left pipette tip racks are placed	[6]
dna_plate_name	Labware containing DNA	"opentrons_96_wellplate_200ul_pcr_full_skirt"
dna_plate_slot	Deck slot of the labware containing DNA	1
cells_plate_name	Labware containing competent cells	"opentrons_96_wellplate_200ul_pcr_full_skirt"
cells_plate_slot	Deck slot of the labware containing competent cells	2
media_plate_name	Labware containing recovery media	"usascientific_12_reservoir_22ml"
media_plate_slot	Deck slot of the labware containing recovery media	3
transformation_plate_name	Labware where transformations occur (located in the thermocycler)	"opentrons_96_wellplate_200ul_pcr_full_skirt"
pre_shock_incubation_temp	Pre heat shock cooling incubation temperature (°C)	4
pre_shock_incubation_time_min	Pre heat shock cooling incubation duration (minutes)	30
heat_shock_temp	Heat shock temperature (°C)	42
heat_shock_time_sec	Heat shock duration (seconds)	30
post_shock_incubation_temp	Post heat shock incubation temperature (°C)	4
post_shock_incubation_time_min	Post heat shock cooling incubation duration (minutes)	2
recovery_temp	Recovery incubation temperature (°C)	37
recovery_time_min	Recovery incubation duration (minutes)	60

Table 1: JSON file configuration parameters for Protocol 1: Heat Shock Transformations (protocol-1-config.json).

Parameter	Description	Example
<code>dna_id</code>	Identifier for the DNA sample	pUC19
<code>dna_source_well</code>	Source well location on the DNA plate containing DNA sample (specify "NA" when adding no DNA)	A1
<code>dna_volume</code>	Volume of DNA to be transferred for transformation (specify 0 when adding no DNA) (μL)	1
<code>cells_id</code>	Identifier for the competent cells	DH5a
<code>cells_source_well</code>	Source well location on the cells plate containing the cells	A1
<code>cells_volume</code>	Volume of competent cells to be transformed (μL)	10
<code>media_id</code>	Identifier for the recovery media	SOC
<code>media_source_well</code>	Source well location on the media plate containing recovery media	A1
<code>media_volume</code>	Volume of media used for recovery (μL)	50
<code>transformation_well</code>	Destination well location on the transformation plate where the transformation mixture (DNA + cells + media) is assembled	A1

Table 2: CSV file configuration data for Protocol 1: Heat Shock Transformations (protocol-1-data.csv).

Example Preparation

Before running the protocol: **1.** Load the P20 8-Channel Pipette (`right_pipette_name`) on the right mount with a compatible Opentrons 96 Tiprack 20 μL (`right_pipette_tiprack_name`) in slot 4 and 5 (`right_pipette_tiprack_slot`). Load the P300 8-Channel Pipette (`left_pipette_name`) on the left mount with a compatible tip rack (`left_pipette_tiprack_name`) in slot 6 (`left_pipette_tiprack_slot`) **2.** Prepare *Opentrons 96 Well* 200 μL plate (`dna_plate_name`) with the plasmid DNA samples and place it on the deck slot 1 (`dna_plate_slot`). **3.** Prepare *Opentrons 96 Well* 200 μL plate (`cells_plate_name`) with the chemically competent *E. coli* cells and keep the cells chilled on ice until needed and then place the plate on the deck slot 2 (`cells_plate_slot`). **4.** Prepare *USA Scientific 12 Well 22 mL* reservoir (`media_plate_name`) with the SOC recovery media and place it on the deck slot 3 (`media_plate_slot`). **5.** Prepare *Opentrons 96 Well* 200 μL plate (`transformation_plate_name`) and place it into the Opentrons thermocycler. **Note:** If using the same labware for multiple reagents, simply repeat the labware name and deck slot in the configuration JSON file.

Example Execution

Step 1: APEX automatically selects which pipettes to use for each transfer, based on available pipettes and required volumes. If there is no pipette in one of the mounts or you don't want APEX to use an installed pipette, specify "NA" for that mount's pipette name parameter.

Step 2: APEX pre-cools the thermocycler to 4 °C (`pre_shock_incubation_temp`). Once the temperature is reached, the thermocycler opens.

Step 3: APEX gently resuspends the cells (`cells_id`) once per source well and distributes the cells

from the cells plate wells A1-H1 (`cells_source_well`) across wells A1-H12 (`transformation_well`) on the transformation plate, with a volume of 10 μL /well (`cells_volume`). APEX will optimise by aspirating sufficient volume for multiple wells when possible and an additional disposal volume equal to the minimum pipette volume is aspirated to ensure accuracy. A single tip is used for all cells transfers, with a blow-out at the source well once finished.

Step 4: APEX transfers 1 μL /well (`dna_volume`) of DNA (`dna_id`) from wells A1-H12 (`dna_source_well`) on the DNA plate into A1-H12 (`transformation_well`) on the transformation plate, mixing gently once. A fresh tip is used for all DNA transfers, and after a blowout, the pipette descends again to add any droplet left at the tip to the well.

Step 5: The thermocycler closes and incubates the cells at 4 °C (`pre_shock_incubation_temp`) for 30 minutes (`pre_shock_incubation_time_min`), followed by a heat shock at 42 °C (`heat_shock_temp`) for 30 seconds (`heat_shock_time_sec`), and an incubation at 4 °C (`post_shock_incubation_temp`) for 2 minutes (`post_shock_incubation_time_min`). APEX automatically calculates the maximum reaction volume based on the content of the well, and configures the thermocycler block accordingly.

Step 6: Thermocycler opens again, and APEX transfers 50 μL /well (`media_volume`) of SOC media (`media_id`) from wells A1-H1 (`media_source_well`) on the media plate into the wells A1-H12 (`destination_well`) on the transformation plate, resuspending the cells 3 times with the mixing volume calculated by APEX as half the total volume of all the contents per well (capped at the pipette's maximum volume).

Step 7: Finally, the thermocycler closes and heats up to 37 °C (`recovery_temp`) and incubates the cells for 60 minutes (`recovery_time_min`). Once finished, the cells are ready for selection via spotting on the selective media plate (APEX Protocol 2)

Protocol 2: Colony Selection

Design

Use the JSON (protocol-2-config.json) and CSV (protocol-2-data.csv) files to create a protocol tailored for colony selection via selective agar plate spotting. Define parameters as described in Table 3 and Table 4 below.

Parameter	Description	Example
pipette_name	Model of the pipette	"p20_multi_gen2"
pipette_mount	Side of the robot the pipette is mounted on	"left" or "right"
tiprack_name	Tip rack type compatible with pipette	"opentrons_96_tiprack_20ul"
tiprack_slot	Deck slot(s) where pipette tip racks are placed	[4]
transformation_plate_name	Labware containing the transformations	"opentrons_96_wellplate_200ul_pcr_full_skirt"
agar_plate_name	Labware of the agar plate	"nunc96grid_96_wellplate_10ul"
agar_plate_area	The base area of the agar plate (cm ²)	9.4692
empty_plate_weight	The weight of the plate without agar (g)	38.92
agar_density	Agar density (g cm ⁻³)	0.911
mixing_volume	Volume of resuspension mixing of the transformed cells before aspirating (μL)	1
additional_volume	Additional aspiration volume ensures accurate dispensing (that volume will not be dispensed onto the plate) (μL)	1
spotting_height	Z-axis offset from agar surface for dispensing (positive values increase height; negative values decrease height; at value 0 the pipette tip will be right at the height of the agar) (mm)	0.5

Table 3: JSON file configuration parameters for Protocol 2: Colony Selection (protocol-2-config.json).

Parameter	Description	Example
sample_id	Identifier for the transformation sample	pUC19-DH5a
agar_plate_slot	Deck slot(s) of the agar plate(s)	1
agar_plate_weight	The weight of the plate(s) with agar (g)	68.92
transformation_source_well	Source well location of the transformed cells	A1
spotting_volume	Volume of cells to spot onto agar (μL)	5
spotting_well	Destination well location on the agar plate where the transformed cells are spotted	A1

Table 4: CSV file configuration data for Protocol 2: Colony Selection (protocol-2-data.csv).

Example Preparation

Before running the protocol: **1.** Load the P20 8-Channel Pipette (pipette_name) on the right mount (pipette_mount) with a compatible Opentrons 96 Tiprack 20 μL (tiprack_name) in slot 4 (tiprack_slot). **2.** Prepare *Opentrons 96 Well* 200 μL plate (transformation_plate_name) with the transformed cells and

place it the Opentrons thermocycler. **3.** Prepare *Nunc OmniTray Single-Well* plate(s) (`agar_plate_name`) with LB agar containing the appropriate selection markers and wait for the agar to solidify and dry. Place it on the deck slot 1 (`agar_plate_slot`). **4.** Record the weights of the plate(s) before (`empty_agar_plate_weight`) and after (`agar_plate_weight`) pouring molten LB agar. **Note:** Proper drying of plates is essential to prevent excessive diffusion of the spotted cells.

Example Execution

Step 1: APEX calculates the agar height based on the plate dimensions (`agar_plate_area`), plate weights (`empty_agar_plate_weight`, `agar_plate_weight`), and agar density (`agar_density`). The user can adjust the spotting height parameter (`spotting_height`) to modify the height of spotting relative to the agar surface. It is recommended to verify the precise plate well depth and agar height using the provided Jupyter notebook.

Step 2: APEX picks up a new tip for each transfer and resuspends the cells 3 times before aspiration with the mixing volume of 20 μL (`mixing_volume`). APEX aspirates 5 μL (`spotting_volume`) of the transformed cells (`sample_id`) from wells A1-H12 (`transformation_source_well`) on the transformation plate and transfers them to the A1-H12 wells (`spotting_well`) on the agar plates. For each transfer an additional volume of 1 μL (`additional_volume`) is aspirated that will not be dispensed. This approach is preferred over using a blowout to prevent sample dispersion and to achieve more precise volume dispensation. After dispensing, the protocol includes a 5-second delay to ensure smooth transfer of the sample.

Step 3: To spot the same sample multiple times without changing tips (technical replicates), multiple destination wells can be listed in the `destination_well` field using the | symbol. For instance, A1 | A2 directs APEX to dispense aliquots from the source well into both wells A1 and A2.

Step 4: Allow the spotted plates to dry before removing them from the OT-2 deck. Transfer the plates for incubation at 37°C for 12-18 hours. After incubation, when colonies have formed, the plates are ready for colony sampling using APEX Protocol 3 (Colony Sampling).

Protocol 3: Colony Sampling

Design

Use the the JSON (protocol-3-config.json) and CSV (protocol-3-data.csv) files to create a protocol tailored for colony sampling. Define parameters as described in Table 5 and Table 6 below.

Parameter	Description	Example
right_pipette_name	Model of the pipette on the right mount (specify "NA" if not used)	"p20_multi_gen2"
right_pipette_tiprack_name	Tip rack type compatible with right pipette	"opentrons_96_tiprack_20ul"
right_pipette_tiprack_slot	Deck slot(s) where right pipette tip racks are placed	[4]
left_pipette_name	Model of the pipette on the left mount (specify "NA" if not used)	"p300_multi_gen2"
left_pipette_tiprack_name	Tip rack type compatible with left pipette	"opentrons_96_tiprack_300ul"
left_pipette_tiprack_slot	Deck slot(s) where left pipette tip racks are placed	[5]
media_plate_name	Labware containing media	"usascientific_12_reservoir_22ml"
media_plate_slot	Deck slot of the labware containing media	1
culture_plate_name	Labware where cells will be propagated	"usascientific_96_wellplate_2.4ml_deep"
culture_plate_slot	Deck slot of the culture plate	2
agar_plate_name	Labware of the agar plate with developed colonies	"nunc96grid_96_wellplate_10ul"
agar_plate_area	The base area of the agar plate (cm ²)	9.4692
empty_plate_weight	The weight of the empty plate without agar (g)	38.92
agar_density	Agar density (g cm ⁻³)	0.911
sampling_method	Method used for sampling colonies from the agar plate	"spiral"
sampling_height	Z-axis offset from agar surface for dispensing (positive values increase height; negative values decrease height; at value 0 the pipette tip will be right at the height of the agar) (mm)	-0.5
spiral_radius	Maximum radius of the spiral sampling pattern (mm)	2.5
spiral_points	Number of discrete sampling points in the spiral path	25
spiral_rotations	Number of complete turns in the spiral sampling pattern	3

Table 5: JSON configuration data details for Protocol 3: Colony Sampling (protocol-3-config.json.)

Example Preparation

Before running the protocol: **1.** Load the P20 8-Channel Pipette (right_pipette_name) on the right mount with a compatible Opentrons 96 Tiprack 20 μ L (right_pipette_tiprack_name) in slot 4 (right_pipette_tiprack_slot). Load the P300 8-Channel Pipette (left_pipette_name) on the left mount with a compatible tip rack (left_pipette_tiprack_name) in slot 5 (left_pipette_tiprack_slot) **2.** Prepare

Parameter	Description	Example
colony_id	Identifier for the bacterial colonies on the agar plate	pUC19-BL21 (DE3)
agar_plate_slot	Deck slot where the agar plate with colonies is placed	3
colony_source_well	Source well location on the agar plate from which to sample colonies	A1
media_id	Identifier of the media	LB-kan
media_source_well	Source well location of the recovery media on the media plate	A3
media_volume	Volume of media for culturing (μL)	500
culture_well	Destination well location on the culture plate where sampled colonies will be inoculated into fresh media	A1

Table 6: CSV file configuration data for Protocol 3: Colony Sampling (protocol-3-data.csv).

USA Scientific 12 Well 22 mL reservoir (media_plate_name) with the LB media and place it on the deck slot 1 (media_plate_slot). **3.** Prepare *USA Scientific 96 Deep Well 2.4 mL* plate (culture_plate_name) and place it on the deck slot 2 (culture_plate_slot). **4.** Prepare *Nunc OmniTray Single-Well* plate(s) (agar_plate_name) with LB agar with developed colonies and place it on the deck slot 3 (agar_plate_slot). Record the weights of the plates(s) with and without LB agar. **5.** Decide on a colony sampling strategy ("spiral" or "pierce").

Example Execution

Step 1: APEX automatically selects which pipettes to use for each transfer, based on available pipettes and required volumes. If there is no pipette in one of the mounts or you don't want APEX to use an installed pipette, specify "NA" for that mount's pipette name parameter.

Step 2: APEX distributes 500 μL /well (media_volume) of growth media supplemented with appropriate antibiotics (media_id) from wells A1-H3 (media_source_well) from the media plate into wells A1-H12 (culture_well) on the culture plate. APEX will optimise by aspirating sufficient volume for multiple wells when possible and an additional disposal volume equal to the minimum pipette volume is aspirated to ensure accuracy. A single tip is used for all media transfers, with a blow-out at the source well once finished.

Step 3: Based on the selected sampling_method, APEX performs colony sampling:

- For "spiral" method: APEX moves the pipette to the center of the specified well containing colonies (colony_id in colony_source_well) and follows a spiral path (defined by spiral_radius, spiral_points, and spiral_rotations) to sample the colonies at sampling_height.
- For "single pierce" method: APEX moves to the center of the specified well and pierces the agar once at sampling_height.

- Step 4:** APEX transfers the sampled cells into the media in the culture plate (culture_well) and mixes the media twice with the mixing volume calculated by APEX as half the total volume of the media volume per well (capped at the pipette's maximum volume).
- Step 5:** After sampling all selected wells with colonies, the culture plate is ready to be moved for shaking incubation at 37°C. After incubation, the plates are ready for protein expression induction using APEX Protocol 4 (Protein Expression Induction).

Protocol 4: Protein Expression Induction

Design

Use the JSON (protocol-4-config.json) and CSV (protocol-4-data.csv) files to create a protocol tailored for protein expression induction. Define parameters as described in Table 7 and Table 8 below.

Parameter	Description	Example
right_pipette_name	Model of the pipette on the right mount (specify "NA" if not used)	"p20_multi_gen2"
right_pipette_tiprack_name	Tip rack type compatible with right pipette	"opentrons_96_tiprack.20ul"
right_pipette_tiprack_slot	Deck slot(s) where right pipette tip racks are placed	[5,6]
left_pipette_name	Model of the pipette on the left mount (specify "NA" if not used)	"p300_multi_gen2"
left_pipette_tiprack_name	Tip rack type compatible with left pipette	"opentrons_96_tiprack.300ul"
left_pipette_tiprack_slot	Deck slot(s) where left pipette tip racks are placed	[7]
media_plate_name	Labware containing media	"usascientific_12_reservoir.22ml"
media_plate_slot	Deck slot of the labware containing media	1
culture_plate_name	Labware where cells will be propagated	"usascientific_96_wellplate.2.4ml_deep"
culture_plate_slot	Deck slot of the culture plate	2
inducer_plate_name	Labware containing the chemical inducer	"opentrons_96_wellplate.200ul_pcr_full_skirt"
inducer_plate_slot	Deck slot of the inducer plate	4
expression_plate_name	Labware where the fresh culture will be propagated and the protein expression will be induced	"usascientific_96_wellplate.2.4ml_deep"
expression_plate_slot	Deck slot of the expression plate	3

Table 7: JSON configuration data details for Protocol 4: Protein Expression Induction (protocol-4-config.json).

Example Preparation

Before running the protocol: **1.** Load the P20 8-Channel Pipette (right_pipette_name) on the right mount with a compatible Opentrons 96 Tiprack 20 μ L (right_pipette_tiprack_name) in slot 5 and 6 (right_pipette_tiprack_slot). Load the P300 8-Channel Pipette (left_pipette_name) on the left mount with a compatible tip rack (left_pipette_tiprack_name) in slot 7 (left_pipette_tiprack_slot) **2.** Prepare *USA Scientific 12 Well 22 mL Reservoir* (media_plate_name) with the media supplemented with antibiotics and place it on the deck slot 1(media_plate_slot). **3.** Prepare *USA Scientific 96 Deep Well 2.4 mL plate* (culture_plate_name) with the cultures and place it on the deck slot 2 (culture_plate_slot). **4.** Prepare *Opentrons 96 Well 200 μ L plate* (inducer_plate_name) with the chemical inducer (e.g. IPTG)

Parameter	Description	Example
culture_id	Identifier for the culture sample	pUC19-BL21 (DE3)
culture_source_well	Source well location on the culture plate containing overnight cultures	A1
culture_volume	Volume of culture to be propagated in fresh media (μL)	5
media_id	Identifier of media	fresh-LB
media_source_well	Source well location on the media plate containing media	A1
media_volume	Volume of media to be distributed per well into the expression plate (μL)	500
inducer_id	Identifier of the chemical inducer	0.1mM-IPTG
inducer_source_well	Source well location on the inducer plate (specify "NA" when adding no inducer)	A1
inducer_volume	Volume of the chemical inducer to be transferred per well into the expression plate (μL) (specify 0 when adding no inducer)	5
expression_well	Destination well location on expression plate where cells will be propagated and protein expression will be induced	A1

Table 8: CSV file configuration data for Protocol 4: Protein Expression Induction (protocol-4-data.csv).

and place it on the deck slot 4 (inducer_plate_slot). **5.** Prepare *USA Scientific 96 Deep Well 2.4 mL* plate (expression_plate_name) as an expression plate where fresh culture will be propagated and protein expression will be induced and place it on the deck slot 3 (expression_plate_slot). **Note:** If using the same labware for multiple reagents, simply repeat the labware name and deck slot in the configuration JSON file.

Example Execution

Step 1: APEX automatically selects which pipettes to use for each transfer, based on available pipettes and required volumes. If there is no pipette in one of the mounts or you don't want APEX to use an installed pipette, specify "NA" for that mount's pipette name parameter.

Step 2: APEX distributes 500 μL/well (media_volume) of LB media supplemented with appropriate antibiotics (media_id) from wells A1-H3 (media_source_well) from the media plate into the wells A1-H12 (expression_well) on the expression plate.

Step 3: APEX aliquots 5 μL (culture_volume) of overnight cultures from wells A1-H12 (culture_source_well) on the culture plate using a new tip each time and inoculates the fresh media in wells A1-H12 (expression_well) on the expression plate by mixing it twice with the mixing volume calculated by APEX as half the total volume of media per well (capped at the pipette's maximum volume).

Step 4: The protocol pauses to allow for incubation. Incubate the plate with shaking at 37°C and

grow culture to desired OD₆₀₀.

Step 5: Upon reaching the desired OD₆₀₀, resume the protocol. APEX adds 5 µL/well (inducer_volume) of IPTG (inducer_id) from wells A1-H1 (inducer_well) on the inducer plate to each culture in wells A1-H12 (expression_well) on the expression plate to induce the protein expression.

Step 6: Once the protocol is finished, continue the incubation of the cultures with shaking at 20 °C to facilitate protein production.

Protocol 5: Colony PCR Sample Preparation

Design

Use the JSON (protocol-5-config.json) and CSV (protocol-5-data.csv) files to create a protocol tailored for colony PCR sample preparation. Define parameters as described in Table 9 and Table 10 below.

Parameter	Description	Example
right_pipette_name	Model of the pipette on the right mount (specify "NA" if not used)	"p20_multi_gen2"
right_pipette_tiprack_name	Tip rack type compatible with right pipette	"opentrons_96_tiprack_20ul"
right_pipette_tiprack_slot	Deck slot(s) where right pipette tip racks are placed	[5]
left_pipette_name	Model of the pipette on the left mount (specify "NA" if not used)	"p300_multi_gen2"
left_pipette_tiprack_name	Tip rack type compatible with left pipette	"opentrons_96_tiprack_300ul"
left_pipette_tiprack_slot	Deck slot(s) where left pipette tip racks are placed	[6]
reagent_plate_name	Labware containing PCR reagents (master mix, primers, water)	"opentrons_96_wellplate_200ul_pcr_full_skirt"
reagent_plate_slot	Deck slot where reagent plate is placed	1
pcr_plate_name	Labware where PCR reactions will be assembled (placed in the thermocycler)	"opentrons_96_wellplate_200ul_pcr_full_skirt"
mmastermix_mix_number	Number of mixing of mastermix with the colony solution	4
amplification_cycles	Number of thermal cycling repetitions	30
lid_temp	Temperature of thermocycler lid (°C)	105
initial_denaturation_temp	Temperature for initial denaturation step (°C)	94
initial_denaturation_time_min	Duration of initial denaturation (minutes)	3
denaturation_temp	Temperature for cycle denaturation step (°C)	94
denaturation_time_sec	Duration of cycle denaturation (seconds)	30
annealing_temp	Temperature for primer annealing step (°C)	55
annealing_time_sec	Duration of primer annealing (seconds)	30
extension_temp	Temperature for DNA extension step (°C)	68
extension_time_sec	Duration of DNA extension (seconds)	60
final_extension_temp	Temperature for final extension step (°C)	68
final_extension_time_min	Duration of final extension (minutes)	5
hold_temp	Final holding temperature (°C)	4

Table 9: JSON file configuration parameters for Protocol 5: Colony PCR Sample Preparation Protocol (protocol-5-config.json).

Example Preparation

Before running the protocol: **1.** Load the P20 8-Channel Pipette (right_pipette_name) on the right mount with a compatible Opentrons 96 Tiprack 20 µL (right_pipette_tiprack_name)

Parameter	Description	Example
colony_id	Identifier for the colony for PCR analysis	sample_1
water_source_well	Source well location in the reagent plate containing nuclease-free water for colony resuspension	A1
water_volume	Volume of nuclease-free water to be dispensed for colony resuspension (μL)	11.5
master_mix_source_well	Source well location in the reagent plate containing the PCR master mix	A2
master_mix_volume	Volume of PCR master mix (μL)	13.5
pcr_well	Destination well location in the PCR plate where the PCR reaction will occur	A1

Table 10: CSV file configuration data for Protocol 5: Colony PCR Sample Preparation Protocol (protocol-5-data.csv).

in slot 5 (right_pipette_tiprack_slot). Load the P300 8-Channel Pipette (left_pipette_name) on the left mount with a compatible tip rack (left_pipette_tiprack_name) in slot 6 (left_pipette_tiprack_slot) **2.** Prepare *Opentrons 96 Well* 200 μL plate (reagent_plate_name) with the PCR master mix and nuclease-free water and place it on the deck slot 2 (reagent_plate_slot). **3.** Prepare *Opentrons 96 Well* 200 μL plate (pcr_plate_name) with where the PCR reaction will occur and place it in the Opentrons thermocycler. Ensure the thermocycler is turned on and available. **4.** Prepare plates with bacterial colonies ready for processing. **5.** Prepare toothpicks or tips for manual colony picking. **Note:** If using the same labware for multiple reagents, simply repeat the labware name and deck slot in the configuration JSON file.

Example Execution

Step 1: APEX automatically selects which pipettes to use for each transfer, based on available pipettes and required volumes. If there is no pipette in one of the mounts or you don't want APEX to use an installed pipette, specify "NA" for that mount's pipette name parameter.

Step 2: APEX distributes 11.5 μL /well (water_volume) of nuclease-free water from wells A1-H1 (water_source_well) on the reagent plate across the wells A1-H12 (pcr_well) on the PCR plate. APEX will optimise by aspirating sufficient volume for multiple wells when possible and an additional disposal volume equal to the minimum pipette volume is aspirated to ensure accuracy. A single tip is used for all cells transfers, with a blow-out at the source well once finished.

Step 3: APEX pauses for manual intervention. The user removes the PCR plate to pick individual bacterial colonies and resuspend them in the distributed water. The user places the plate into the thermocycler module and resumes the protocol.

Step 4: APEX transfers 13.5 μ L (`mastermix_volume`) of the PCR master mix from the well A2-H2 (`mastermix_source_well`) on the reagent plate into each well containing resuspended colonies (`pcr_well`) on the PCR plate. The contents of each well are mixed 4 times (`mastermix_mix_number`) with the mixing volume calculated by APEX as half the total volume of master mix and water per well (capped at the pipette's maximum volume). A new tip is used for each transfer with blow-out at the destination well.

Step 5: The thermocycler closes its lid and sets lid temperature to 105 °C (`lid_temp`). Initial denaturation is performed at 94 °C (`initial_denaturation_temp`) for 3 minutes (`initial_denaturation_time_min`). The PCR thermal cycling then begins, with 30 repetitions (`amplification_cycles`) of: denaturation at 94 °C (`denaturation_temp`) for 30 seconds (`denaturation_time_sec`), primer annealing at 55 °C (`annealing_temp`) for 30 seconds (`annealing_time_sec`), and DNA extension at 68 °C (`extension_temp`) for 60 seconds (`extension_time_sec`). APEX automatically calculates the maximum reaction volume based on the combined PCR master mix, and colony suspension in water volumes, and configures the thermocycler block accordingly. After cycling completes, a final extension at 68 °C (`final_extension_temp`) for 5 minutes (`final_extension_time_min`). Samples are then held at 4 °C (`hold_temp`) until retrieval and are ready for subsequent analysis via gel electrophoresis.

Protocol 6: SDS PAGE Sample Preparation

Design

Use the JSON (protocol-6-config.json) and CSV (protocol-6-config.csv) files to create a protocol tailored for SDS PAGE sample preparation. Define parameters as described in Table 11 and Table 12 below.

Parameter	Description	Example
right_pipette_name	Model of the pipette on the right mount (specify "NA" if not used)	"p20_multi_gen2"
right_pipette_tiprack_name	Tip rack type compatible with right pipette	"opentrons_96_tiprack_20ul"
right_pipette_tiprack_slot	Deck slot(s) where right pipette tip racks are placed	[5]
left_pipette_name	Model of the pipette on the left mount (specify "NA" if not used)	"p300_multi_gen2"
left_pipette_tiprack_name	Tip rack type compatible with left pipette	"opentrons_96_tiprack_300ul"
left_pipette_tiprack_slot	Deck slot(s) where left pipette tip racks are placed	[6]
culture_plate_name	Labware containing culture	"usascientific_96_wellplate_2.4ml_deep"
culture_plate_slot	Deck slot where culture plate is placed	1
reagent_plate_name	Labware containing reagents (lysis buffer, SDS buffer)	"usascientific_12_reservoir_22ml"
reagent_plate_slot	Deck slot where reagent plate is placed	2
lysate_plate_name	Labware to which lysate will be transferred	"usascientific_12_reservoir_22ml"
lysate_plate_slot	Deck slot where lysate plate is placed	3
destination_plate_name	Labware used for SDS sample prep (located in the thermocycler)	"opentrons_96_wellplate_200ul_pcr_full_skirt"
culture_mix_number	Number of mixing repetitions before culture transfer	3
culture_mix_volume	Volume for mixing culture (μL)	20

Table 11: JSON file configuration parameters for Protocol 6: SDS PAGE Sample Preparation (protocol-6-config.json).

Parameter	Description	Example
sample_id	Identifier for the sample	sample_1
culture_source_well	Source well location in the culture plate containing bacterial culture	A1
culture_volume	Volume of bacterial culture to be analysed for expression (μL)	150
lysis_buffer_source_well	Source well location in the reagent plate containing lysis buffer	A1
lysis_buffer_volume	Volume of lysis buffer (μL)	40
sds_buffer_source_well	Source well location in the reagent plate containing SDS buffer	A1
sds_buffer_volume	Volume of SDS buffer (μL)	10
destination_well	Destination well location in the destination plate where the SDS sample prep will occur	A1

Table 12: CSV file configuration data for Protocol 6: SDS PAGE Sample Preparation (protocol-6-data.csv).

Example Preparation

Before running the protocol: **1.** Load the P20 8-Channel Pipette (`right_pipette_name`) on the right mount with a compatible Opentrons 96 Tiprack 20 μ L (`right_pipette_tiprack_name`) in slot 5 (`right_pipette_tiprack_slot`). Load the P300 8-Channel Pipette (`left_pipette_name`) on the left mount with a compatible tip rack (`left_pipette_tiprack_name`) in slot 6 (`left_pipette_tiprack_slot`) **2.** Prepare *USA Scientific 96 Deep Well 2.4 mL* plate (`culture_plate_name`) with the bacterial cultures and place it on the deck slot 1 (`culture_plate_slot`). **3.** Prepare *USA Scientific 12 Well 22 mL* reservoir (`reagent_plate_name`) with the reagents (lysis buffer and SDS buffer) and place it on the deck slot 2 (`reagent_plate_slot`). **4.** Prepare *USA Scientific 96 Deep Well 2.4 mL* plate (`lysate_plate_name`) where lysate will be transferred and place it on the deck slot 3 (`culture_plate_slot`). **5.** Prepare *Opentrons 96 Well 200 μ L* plate (`destination_plate_name`) where the SDS sample prep will occur and place it on the int he Opentrons thermocycler. Ensure the thermocycler is turned on and available. **6.** Set up the external centrifuge with plate rotors. **Note:** If using the same labware for multiple reagents, simply repeat the labware name and deck slot in the configuration JSON file.

Example Execution

Step 1: APEX automatically selects which pipettes to use for each transfer, based on available pipettes and required volumes. If there is no pipette in one of the mounts or you don't want APEX to use an installed pipette, specify "NA" for that mount's pipette name parameter.

Step 2: APEX transfers 150 μ L (`culture_volume`) of bacterial cultures from the wells A1-H12 (`culture_source_well`) on the culture culture plate to the wells A1-H12 (`destination_well`) on the lysis plate. Before aspiration, APEX mixes the cultures 3 times (`culture_mix_number`) with a volume of 100 μ L (`culture_mix_volume`) to resuspend the cells. A new tip is used for each transfer with blow-out at the destination well.

Step 3: The protocol pauses with a message prompting the user for manual centrifugation of the lysis plate. After centrifugation and returning the plate to the deck, protocol is resumed and pipetting parameters are automatically adjusted to avoid disturbing the cell pellet. The pipette's well bottom clearance is set to 0.5 mm (`pellet_well_bottom_clearance`) which can be adjusted based on the pellet size. Aspiration flow rate is reduced to 25 % (`supernatant_flow_rate_aspirate`) for gentle supernatant removal. The supernatant is discarded into the waste container, and pipette settings return to default values.

- Step 4:** APEX transfers 40 μL /well (lysis_buffer_volume) of the lysis buffer from wells A1-H1 (lysis_buffer_source_well) on the reagent plate to each well containing the cell pellets (destination_well). APEX performs 5 mixes (lysis_mix_number) to thoroughly resuspend the pellet. Mixing volume is calculated by APEX as half the lysis buffer volume per well (capped at the maximum volume of the pipette). A new tip is used for each transfer with blow-out at the destination well.
- Step 5:** Samples are incubated for 30 min (lysis_incubation_time_min) at room temperature. The protocol then pauses for manual centrifugation of the lysis plate to pellet cell debris.
- Step 6:** APEX transfers 40 μL (lysis_buffer_volume) of the lysate from wells A1-H12 (destination_well) on the lysate plate into the wells A1-H12 (destination_well) on the destination plate located in the thermocycler. A new tip is used for each transfer with blow-out at the destination well.
- Step 7:** APEX transfers 10 μL /well (sds_buffer_volume) of the 5x SDS buffer from wells A2-H12 (sds_buffer_source_well) on the reagent plate to each well A1-H12 (destination_well) containing the lysate. APEX performs 4 mixes (sds_mix_number). Mixing volume is calculated as half the total volume of lysate and SDS buffer in the well (capped at the maximum volume of the pipette). A new tip is used for each transfer with a blow-out at the destination well.
- Step 8:** The thermocycler closes its lid and sets the lid temperature to 105 °C (lid_temp). The maximum reaction volume is calculated by APEX and thermocycler block max volume is set accordingly. Samples are heated to 95 °C (denaturation_temp) for 5 min (denaturation_time_min). The thermocycler deactivates both lid and block heater once finished. The denatured protein samples are now ready for gel loading.

Protocol 7: OD₆₀₀ Reading Sample Preparation

Design

Use the JSON (protocol-6-config.json) and CSV (protocol-7-config.csv) files to create a protocol tailored for OD₆₀₀ reading sample preparation. Define parameters as described in 13 and 14 below.

Parameter	Description	Example
right_pipette_name	Model of the pipette on the right mount (specify "NA" if not used)	"p20_multi_gen2"
right_pipette_tiprack_name	Tip rack type compatible with right pipette	"opentrons.96_tiprack_20ul"
right_pipette_tiprack_slot	Deck slot(s) where right pipette tip racks are placed	[5]
left_pipette_name	Model of the pipette on the left mount (specify "NA" if not used)	"p300_multi_gen2"
left_pipette_tiprack_name	Tip rack type compatible with left pipette	"opentrons.96_tiprack_300ul"
left_pipette_tiprack_slot	Deck slot(s) where left pipette tip racks are placed	[6]
diluent_plate_name	Labware containing diluent (e.g. water, media)	"corning_96_wellplate_360ul_flat"
diluent_plate_slot	Deck slot where diluent plate is placed	1
culture_plate_name	Labware containing culture	"usascientific_96_wellplate_2.4ml_deep"
culture_plate_slot	Deck slot where culture plate is placed	2
reading_plate_name	Labware for reading OD ₆₀₀ of the cultures	"opentrons.96_wellplate_200ul_pcr_full_skirt"
reading_plate_slot	Deck slot where reading plate is placed	3
culture_mix_number	Number of mixing of bacterial culture before transfer	3
culture_mix_volume	Volume for mixing bacterial culture (μL)	20

Table 13: JSON file configuration parameters for Protocol 7: OD₆₀₀ Reading Sample Preparation (protocol-7-config.json).

Parameter	Description	Example
culture_id	Identifier for the sample	sample_1
diluent_source_well	Source well location in the diluent plate containing diluent	A1
diluent_volume	Volume of diluent to be dispensed for diluting cultures (μL)	90
culture_source_well	Source well location in the culture plate containing bacterial culture (specify "NA" when adding no culture for the blank)	A1
culture_volume	Volume of bacterial culture to be transferred (specify 0 when adding no culture for the blank) (μL)	10
reading_well	Destination well position in the reading plate where the bacterial culture will be diluted for OD ₆₀₀ measurement	A1

Table 14: CSV file configuration data for Protocol 7: OD₆₀₀ Reading Sample Preparation (protocol-7-data.csv).

Example Preparation

Before running the protocol: **1.** Load the P20 8-Channel Pipette (`right_pipette_name`) on the right mount with a compatible Opentrons 96 Tiprack 20 μL (`right_pipette_tiprack_name`) in slot 5 (`right_pipette_tiprack_slot`). Load the P300 8-Channel Pipette (`left_pipette_name`) on the left mount with a compatible tip rack (`left_pipette_tiprack_name`) in slot 6 (`left_pipette_tiprack_slot`) **2.** Prepare *USA Scientific 12 Well 2.4 mL* plate (`culture_plate_name`) with the cultures and place it on the deck slot q (`culture_plate_slot`). **3.** Prepare *Corning 96 Well 360 μL* plate (`reading_plate_name`) and place it on the deck slot 3 (`reading_plate_slot`). **4.** Prepare USA Scientific 12 Well 22 mL Reservoir (`diluent_plate_name`) with the diluent (typically water or media) and place it on the deck slot 2 (`diluent_plate_slot`). **Note:** If using the same labware for multiple reagents, simply repeat the labware name and deck slot in the configuration JSON file.

Example Execution

Step 1: APEX automatically selects which pipettes to use for each transfer, based on available pipettes and required volumes. If there is no pipette in one of the mounts or you don't want APEX to use an installed pipette, specify "NA" for that mount's pipette name parameter.

Step 2: APEX distributes 90 μL /well (`diluent_volume`) of diluent from wells A1-H1 (`diluent_source_well`) on the diluent plate across the wells A1-H11 (`reading_well`) and 100 μL /well (`diluent_volume` + `culture_volume`) in wells A12-H12. APEX will optimise by aspirating sufficient volume for multiple wells when possible and an additional disposal volume equal to the minimum pipette volume is aspirated to ensure accuracy. A single tip is used for all diluent transfers, with a blow-out at the source well once finished.

Step 3: APEX transfers 10 μL (`culture_volume`) of bacterial cultures from wells A1-H12 (`culture_source_well`) on the culture plate to the wells A1-H11 (`reading_well`) on the reading plate. Before aspiration, APEX mixes the cultures 5 times (`culture_mix_number`) with a volume of 250 μL (`culture_mix_volume`) to resuspend the cells. A new tip is used for each transfer with a blow-out at the destination well.

Step 4: For calibration blanks, specify "NA" in the `culture_source_well` field and "0" in the `culture_source_volume`. For blank samples, set `diluent_volume` equal to the total desired final volume. Once the protocol is finished, samples can be analysed using a plate reader to determine bacterial cell density.